



Alere Triage® Solutions Product Catalog

Alere Triage® Solutions Product List

Product Name	Alere PN	Quantity & Format
Alere Triage® BNP Test		
Alere Triage® BNP Test	98010XR	15 Tests
Alere Triage® BNP Test	98000XR	25 Tests
Alere Triage® BNP Control 1	98013XR	5 x 0.25mL Vials
Alere Triage® BNP Control 2	98014XR	5 x 0.25mL Vials
Alere Triage® BNP Calibration Verification	98015XR	2 x 5 Levels
Alere Triage® BNP Test for the Beckman Coulter Access Family of Immunoassay Systems		
Alere Triage® BNP Reagent Packs	98200	2 x 50 Tests
Alere Triage® BNP QC Controls	98201	2 x 3 Levels
Alere Triage® BNP Calibrators	98202	1 x 6 Levels
Alere Triage® D-Dimer Test		
Alere Triage® D-Dimer Test	98100	25 Tests
Alere Triage® Total 5 Control 1	88753	5 x 0.25mL Vials
Alere Triage® Total 5 Control 2	88754	5 x 0.25mL Vials
Alere Triage® Total 5 Calibration Verification	88755	2 x 5 Levels
Alere Triage® Cardiac Panel		
Alere Triage® Cardiac Panel	97000HS	25 Tests
Alere Triage® Cardiac Panel (CK-MB and Troponin I)	97022HS	25 Tests
Alere Triage® Cardiac Panel (Troponin I)	97021HS	25 Tests
Alere Triage® Total 5 Control 1	88753	5 x 0.25mL Vials
Alere Triage® Total 5 Control 2	88754	5 x 0.25mL Vials
Alere Triage® Total 5 Calibration Verification	88755	2 x 5 Levels
Alere Triage® TOX Drug Screen		
Alere Triage® TOX Drug Screen	94400	25 Tests
Alere Triage® TOX Drug Screen-10	94401	25 Tests
Alere Triage® TOX Drug Screen-9	94402	25 Tests
Alere Triage® TOX Drug Screen Control 1	94413	5 x 0.25mL Vials
Alere Triage® TOX Drug Screen Control 2	94414	5 x 0.25mL Vials
Alere Triage® Meter		
Alere Triage® MeterPro	55070	1
Alere Triage® MeterPro User Manual	26213en	1
Alere Triage® Supervisor CODE CHIP™ Module	80035A	1
Alere Triage® Laser Barcode Scanner	52111	1
Alere Triage® Meter Printer Thermal Paper Roll	52031	1
Alere Triage® Transfer Pipettes	52040	25/Bag
Alere Triage® Data Management/Connectivity		
Alere Triage Census® ATM Software Kit (includes all the items listed below)	55069	1
Alere Triage Census® ATM Software Version 4.0 User's Guide	22643	1
Alere Triage Census® ATM Software Version 4.0 CD	52140	1
Alere Triage® Meter Serial Interface Cable (RS-232 DB9M-DB9F)	52039	1

Tests Available on the Alere Triage® Meter

Test	Intended Use
Alere Triage® BNP Test	The Alere Triage® BNP Test is a rapid, point of care fluorescence immunoassay to be used with the Alere Triage® Meters for the quantitative measurement of B-type natriuretic peptide (BNP) in EDTA anticoagulated whole blood or plasma specimens. The test is intended to be used as an aid in the diagnosis and assessment of severity of congestive heart failure (also referred to as heart failure). The test also is used for the risk stratification of patients with acute coronary syndromes and for the risk stratification of patients with heart failure.
Alere Triage® D-Dimer Test	The Alere Triage® D-Dimer Test is a fluorescence immunoassay to be used with the Alere Triage® Meters for the quantitative determination of cross-linked fibrin degradation products containing D-dimer in EDTA anticoagulated whole blood and plasma specimens. The test is used as an aid in the assessment and evaluation of patients suspected of having disseminated intravascular coagulation or thromboembolic events including pulmonary embolism.
Alere Triage® Cardiac Panel	The Alere Triage® Cardiac Panel is a fluorescence immunoassay to be used with the Alere Triage® Meters for the quantitative determination of creatine kinase MB (CK-MB), myoglobin and troponin I in EDTA anticoagulated whole blood or plasma specimens. The test is used as an aid in the diagnosis of myocardial infarction (injury).
Alere Triage® TOX Drug Screen	The Alere Triage® TOX Drug Screen is a fluorescence immunoassay intended to be used with the Alere Triage® Meters for the qualitative determination of the presence of drug and/or the major metabolites above the threshold concentrations of up to 10 distinct drug classes, including assays for acetaminophen/paracetamol, amphetamines, methamphetamines, barbiturates, benzodiazepines, cocaine, methadone, opiates, phencyclidine, THC and tricyclic antidepressants in urine. The acetaminophen/paracetamol assay will yield positive results when acetaminophen/paracetamol is ingested at or above therapeutic doses*.

*Please refer to the Alere Triage® Package Insert for the complete intended use statement.



The Alere Triage® MeterPro combines sophisticated technology with user-friendly ergonomics to deliver high-performance testing in the laboratory or at the point of care (POC).

- Results in about 15-20 minutes
- Whole blood, plasma, or urine
- Largest immunoassay menu available in a POC platform
- Fluorescence spectrophotometer analyzer
- Customizable test panels via Test Select™ feature

Alere Triage® Product Analytes

	Myoglobin	CK-MB	Troponin I	BNP	D-dimer	TOX
Alere Triage® Cardiac Panel	■	■	■			
Alere Triage® Cardiac Panel (CK-MB and Troponin I)		■	■			
Alere Triage® Cardiac Panel (Troponin I)			■			
Alere Triage® BNP Test				■		
Alere Triage® D-Dimer Test					■	
Alere Triage® TOX Drug Screen APAP, AMP, mAMP, BAR, BZO, COC, MTD, OPI, PCP, THC, TCA						■
Alere Triage® TOX Drug Screen-10 AMP, mAMP, BAR, BZO, COC, MTD, OPI, PCP, THC, TCA						■
Alere Triage® TOX Drug Screen-9 AMP, mAMP, BAR, BZO, COC, OPI, PCP, THC, TCA						■



Call **1.877.441.7440**
or visit **alere.com** today