

COUNTER CO

- **CO Check pro**, analyzer suitable for measurement of carbon monoxide, ideal for use in clinics for smoking cessation.

The results are exported to % COHb (percent concentration of carbon monoxide in hemoglobin) or PPM (Parts per million) and has additional LED lighting according to their coloring, warn the user, Non-smoker, Light, Medium or Heavy smoker.

Measuring time : 15 sec

Required the reservation of respiration for 15 sec and after full exhalation 10 sec

CONSUMABLE

MATERIALS

Paper valves one-way valves

The meter automatically compensates for the level of CO in the atmosphere.



To keep the meter in excellent condition, require six monthly calibration through appropriate mixture gas